



**瓶口分液器**  
**Bottle-top Dispenser**

# Bottle-top Dispenser

## 产品特点

- 具有极强的化学耐受性
- 可高温灭菌消毒
- 可调节控制旋钮使量程调节更为便捷
- 采用PTFE FEP BSG PP 等材料制成
- 最大耐压500mbar, 最大耐粘性500mm<sup>2</sup>/s, 最大耐液体温度为40, 最大耐液体密度 2.2g/cm<sup>2</sup>

## Product Features

- Excellent chemical resistance
- Fully autoclavable at 121°C
- Adjustable control knob makes setting the volume easy
- Made of PTFE, FEP, BSG, PP
- Vapor pressure Max. 500mbar, viscosity max. 500mm<sup>2</sup> / s, temperature max. 40°C, density max. 2.2g/cm<sup>2</sup>

## 订购详情 Ordering Information

产品编号	量程	最小分度	准确度	精确度
Product Code	Volume Range	Subdivision	Accuracy	Precision
30003010	1 mL	—	±2.00% 20ul	1.00% 10ul
30003020	0.5-5 mL	0.1 mL	±0.50% 25ul	0.10% 5ul
30003030	1.0-10 mL	0.2 mL	±0.50% 50ul	0.10% 10ul
30003040	2.5-25 mL	0.5 mL	±0.50% 125ul	0.10% 25ul
30003050	5.0-50 mL	1 mL	±0.50% 250ul	0.10% 50ul

## 附件 Accessories

产品编号	描述	产品编号	描述
Product Code	Describe	Product Code	Describe
30003110	瓶口适配器 Adapt 45/40mm	30003150	瓶口适配器 Adapt 32/28mm
30003120	瓶口适配器 Adapt 45/32mm	30003160	棕色试剂瓶(1L) Bottle (Brown, 1L)
30003130	瓶口适配器 Adapt 45/38mm	30003170	进液管(米) Filling tube (m)
30003140	瓶口适配器 Adapt 32/25mm	30003180	出液管(米) Discharge tube (m)



**RONGTAI**  
BIOCHEMICAL

上海荣泰生化工程有限公司  
SHANGHAI RONGTAI BIOCHEMICAL ENGINEERING CO., LTD.

Email: [sales@rt-bio.com](mailto:sales@rt-bio.com) Website: [www.rt-bio.com](http://www.rt-bio.com)

地址: 上海市松江区曹农路255号8幢  
邮编: 201612  
电话: 8621 3701 3386 • 3701 3385  
传真: 8621 5988 2370 企业QQ: 800033773

Address: Building 8, No 255, Cao Nong Road,  
Song Jiang District, Shanghai, 201612, P. R. China  
Tel: 8621 3701 3386 • 3701 3385  
Fax: 8621 5988 2370 QQ EIM: 800033773