

RONGTAI
BIOCHEMICAL



上海荣泰生化工程有限公司
SHANGHAI RONGTAI BIOCHEMICAL ENGINEERING CO., LTD.



ePipette

ePipette电动移液器

- 电动移液器，双屏设计，可充电锂电池，量程范围0.1-10000ul。
- 超轻设计，流线造型，符合人体工程学。
- 大液晶屏位于枪身部分，锁定感应设计，USB充电模式。
- 组合旋钮设计，所有操作均可通过组合旋钮来完成，只需一个按键便可完成全部吸液和排液，设置容量等操作。
- 吸液和排液提供9个档位的调节速度。
- 精密度和准确度高，保证移液的准确性和平行性。

Key Features

High accuracy, high performance stepper motor ensures accuracy and repeatability, eliminating manual pipetting errors

- 2 LCM Display, Covered 0.1-10000ul volume range.
- A motor driven digital control pipette with multifunction.
- A special locking induction design, 2 buttons handle all operational settings.
- Light weight, ergonomic design, small body dimension for easy handling that guarantees fatigue free pipetting
- Pipetting, Mixing, Stepper and Dilution.
- Adjustable speeds for aspiration and dispensing.
- Li-ion battery and dual charging modes enable longer operation time.

订购详情 Ordering Information

产品编号	品名	量程 (ul)	增量	测试点	准确度	精确度	适用管嘴代码
Product Code	Description	Volume Range	Increment	Test Volume	Accuracy	Precision	Tips Code
90201010	ePipette	0.1-10ul	0.01ul	10ul	± 0.90%	≤ 0.40%	80102010
				5ul	± 1.00%	≤ 0.70%	
				1ul	± 2.50%	≤ 1.50%	
90201020	ePipette	5-120ul	0.10ul	120ul	± 0.40%	≤ 0.15%	80102060
				60ul	± 0.60%	≤ 0.20%	
				12ul	± 2.00%	≤ 1.00%	
90201030	ePipette	10-300ul	0.20ul	300ul	± 0.40%	≤ 0.15%	80102060
				150ul	± 0.60%	≤ 0.20%	
				30ul	± 1.50%	≤ 0.80%	
90201040	ePipette	50-1000ul	1.00ul	1000ul	± 0.40%	≤ 0.15%	80102090
				500ul	± 0.60%	≤ 0.20%	
				100ul	± 1.50%	≤ 0.50%	
90201050	ePipette	100-5000ul	5.00ul	5000ul	± 0.50%	≤ 0.15%	80102100
				2500ul	± 0.80%	≤ 0.20%	
				500ul	± 1.00%	≤ 0.40%	
90201060	ePipette	500-10000ul	10.00ul	10000ul	± 0.60%	≤ 0.20%	80102153
				5000ul	± 1.20%	≤ 0.30%	
				1000ul	± 3.00%	≤ 0.60%	

RONGTAI
BIOCHEMICAL

上海荣泰生化工程有限公司
SHANGHAI RONGTAI BIOCHEMICAL ENGINEERING CO., LTD.

Email: sales@rt-bio.com Website: www.rt-bio.com

地址: 上海市松江区曹农路255号8幢
邮编: 201612
电话: 8621 3701 3386 • 3701 3385
传真: 8621 5988 2370 企业QQ: 800033773

Address: Building 8, No 255, Cao Nong Road,
Song Jiang District, Shanghai, 201612, P. R. China
Tel: 8621 3701 3386 • 3701 3385
Fax: 8621 5988 2370 QQ EIM: 800033773