



多道移液器
Multi-channel Pipette

PIPETTE

Multi-channel Pipette

手动多道可调移液器的特点:

- 管嘴推出器可同时推出多道吸嘴，高效省力
- 各量程的多道移液器均适用于标准96孔板
- 液头可360度旋转，方便移液
- 每道管嘴连件都有独立的活塞装置，维修保养便捷
- 特别的管嘴连件设计，易于观察吸嘴的密封状况

Features for Multi-channel Pipette

- 8 channel pipettes are available for standard 96-well plate
- Dispensing head rotates for optimum pipetting convenience
- Individual piston and tip cone assemblies allow easy repair and maintenance
- Compound material tip cone design allows visual seal verification
- Can be used with universal style pipette tips

订购详情 Ordering Information

产品编号	通道数	量程	增量	测试点	准确度	精确度	适用管嘴代码
Product Code	Channels	Volume Range	Increment	Test Volume	Accuracy	Precision	Tips Code
20001110	8-channel	0.5-10ul	0.1ul	10ul	1.50%	1.50%	80102010
				5ul	2.50%	2.50%	
				1ul	4.00%	4.00%	
20001115	8-channel	2-20ul	0.5ul	20ul	0.90%	0.40%	80102060
				10ul	1.20%	1.00%	
				2ul	3.00%	2.00%	
20001120	8-channel	5-50ul	0.5ul	50ul	1.00%	0.50%	80102060
				25ul	1.50%	1.00%	
				5ul	3.00%	2.00%	
20001121	8-channel	10-100ul	1ul	100ul	0.80%	0.15%	80102060
				50ul	1.00%	0.40%	
				10ul	3.00%	1.50%	
20001122	8-channel	20-200ul	1ul	200ul	0.60%	0.15%	80102060
				100ul	0.80%	0.30%	
				20ul	3.00%	1.00%	
20001129	8-channel	30-300ul	5ul	300ul	0.70%	0.25%	80102060
				150ul	1.00%	0.50%	
				30ul	1.50%	0.80%	
20001130	8-channel	50-300ul	5ul	300ul	0.70%	0.25%	80102060
				150ul	1.00%	0.50%	
				50ul	1.50%	0.80%	
20001210	12-channel	0.5-10ul	0.1ul	10ul	1.50%	1.50%	80102010
				5ul	2.50%	2.50%	
				1ul	4.00%	4.00%	
20001215	12-channel	2-20ul	0.5ul	20ul	0.90%	0.40%	80102060
				10ul	1.20%	1.00%	
				2ul	3.00%	2.00%	
20001220	12-channel	5-50ul	0.5ul	50ul	1.00%	0.50%	80102060
				25ul	1.50%	1.00%	
				5ul	3.00%	2.00%	
20001221	12-channel	10-100ul	1ul	100ul	0.80%	0.15%	80102060
				50ul	1.00%	0.40%	
				10ul	3.00%	1.50%	
20001222	12-channel	20-200ul	1ul	200ul	0.60%	0.15%	80102060
				100ul	0.80%	0.30%	
				20ul	3.00%	1.00%	
20001229	12-channel	30-300ul	5ul	300ul	0.70%	0.25%	80102060
				150ul	1.00%	0.50%	
				30ul	1.50%	0.80%	
20001230	12-channel	50-300ul	5ul	300ul	0.70%	0.25%	80102060
				150ul	1.00%	0.50%	
				50ul	1.50%	0.80%	

RONGTAI
BIOCHEMICAL

上海荣泰生化工程有限公司
SHANGHAI RONGTAI BIOCHEMICAL ENGINEERING CO., LTD.

Email: sales@rt-bio.com Website: www.rt-bio.com

地址: 上海市松江区曹农路255号8幢
邮编: 201612
电话: 8621 3701 3385 ▪ 8621 3701 3386
传真: 8621 5988 2370 企业QQ: 800033773

Address: Building 8, No 255, Cao Nong Road,
Song Jiang District, Shanghai, 201612, P. R. China
Tel: 8621 3701 3385 ▪ 8621 3701 3386
Fax: 8621 5988 2370 QQ EIM: 800033773