



单道移液器
Single Channel Digital Variable Pipette

PIPETTE

Single Channel Digital Variable Pipette

半只消毒手动单道可调式移液器的特点：

- 人体学的手柄挂钩设计，可轻松挂靠在手掌上
- 管嘴弹出按钮设计轻巧
- 清晰的显示读数窗口
- 管嘴推出器可拆卸，可高温高压消毒，便于清洗和保养
- 备有配套的移液器挂钩、拆卸工具

PIPETTE

Features for Single Channel Digital Variable Pipette

- An ergonomic design, convenience of use and simplicity of operation
- Ease of calibration and maintenance
- Availability for autoclave of part of the device
- Detailed instruction in Chinese and English
- A hook, disassembly tool and consumable replacement parts

订购详情 Ordering Information

产品编号	通道数	量程	增量	测试点	准确度	精确度	适用管嘴代码
Product Code	Channels	Volume Range	Increment	Test Volume	Accuracy	Precision	Tips Code
20001005	1	0.1-2.5ul	0.05ul	2.5ul 1.25ul 0.25ul	2.50% 3.00% 12.00%	2.00% 3.00% 6.00%	80102010
20001010	1	0.5-10ul	0.1ul	10ul 5ul 1ul	1.00% 1.50% 2.50%	0.80% 1.50% 1.50%	80102010
20001015	1	2-20ul	0.5ul	20ul 10ul 2ul	0.90% 1.20% 3.00%	0.40% 1.00% 2.00%	80102060
20001020	1	5-50ul	0.5ul	50ul 25ul 5ul	0.60% 0.90% 2.00%	0.30% 0.60% 2.00%	80102060
20001025	1	10-100ul	1ul	100ul 50ul 10ul	0.80% 1.00% 3.00%	0.15% 0.40% 1.50%	80102010
20001030	1	20-200ul	1ul	200ul 100ul 20ul	0.60% 0.80% 3.00%	0.15% 0.30% 1.00%	80102060
20001040	1	100-1000ul	5ul	1000ul 500ul 100ul	0.60% 0.70% 2.00%	0.20% 0.25% 0.70%	80102090
20001050	1	1-5ml	0.05ml	5ml 2.5ml 1ml	0.50% 0.60% 0.70%	0.15% 0.30% 0.30%	80102100
20001060	1	1-10ml	0.1ml	10ml 5ml 1ml	0.60% 1.20% 3.00%	0.20% 0.30% 0.60%	80102135

RONGTAI
BIOCHEMICAL

上海荣泰生化工程有限公司
SHANGHAI RONGTAI BIOCHEMICAL ENGINEERING CO., LTD.

Email: sales@rt-bio.com Website: www.rt-bio.com

地址：上海市松江区曹农路255号8幢
邮编：201612
电话：8621 3701 3386 ▪ 3701 3385
传真：8621 5988 2370 企业QQ：800033773

Address: Building 8, No 255, Cao Nong Road,
Song Jiang District, Shanghai, 201612, P. R. China
Tel: 8621 3701 3386 ▪ 3701 3385
Fax: 8621 5988 2370 QQ EIM: 800033773